

RM3479 Germanium metal granules**Product Number Packing**

RM3479 : 1G

Product Information

Product Code : RM3479
Product Name : Germanium metal granules
Atomic Formula : Ge
Atomic Weight : 72.64
CAS No. : 7440-56-4
EC No. : 231-164-3
HS Code : 8112 92 95
Shelf Life : 4 years

Technical Specification

Appearance : Metallic grey granules or chips
Solubility : Insoluble in water
Particle size : 0.7 - 3.5 mm
ICP : Confirms Germanium component
Trace Metal Analysis : <= 20.0 ppm
Arsenic (As) : ppm
Boron (B) : ppm
Copper (Cu) : ppm
Cobalt (Co) : ppm
Lead (Pb) : ppm
Lithium (Li) : ppm
Manganese (Mn) : ppm
Molybdenum (Mo) : ppm
Nickel (Ni) : ppm
Zinc (Zn) : ppm
Assay (Based on Trace Metal Analysis) : min. 99.99 %

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37/38
S-Phrase(s) : 26-36/39

WGK : 3
RTECS : LY5200000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods