

RM5366 DL-Noradrenaline**Product Number Packing**

RM5366 : 1G

Product Information

Product Code : RM5366
Product Name : DL-Noradrenaline
Synonym : DL-Norepinephrine
Molecular Formula : $C_8H_{11}NO_3$
Molecular Weight : 169.18
CAS No. : 138-65-8
EC No. : 205-337-9
Shelf Life : 4 years

Technical Specification

Appearance : White to brown crystals or powder or solid
Solubility : 10 mg soluble in 10 mL of water
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 185 - 195°C (dec.)
Assay (NT) : 98.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H302-H318
Precautionary Statement(s) : P280-P305+P351+P338
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 28-41
S-Phrase(s) : 28-36/37/39-45
WGK : 3
Storage Temperature(°C) : Store at 2 - 8°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods