

## RM6249 Agarose, medium EEO For electrophoresis

### Product Number Packing

---

RM6249	: 100G
RM6249	: 10G
RM6249	: 1KG
RM6249	: 25G

### Product Information

---

Product Code	: RM6249
Product Name	: Agarose, medium EEO For electrophoresis
CAS No.	: 9012-36-6
Shelf Life	: 4 years

### Technical Specification

---

Appearance	: White crystalline powder.
Solubility	: Soluble in hot, boiling water.
Melting point	: $87.5 \pm 2.5^{\circ}\text{C}$
Gelling strength (1.0%)	: $> 1000 \text{ g/ cm}^2$
Gel point (1.5%)	: $33 \pm 1.5^{\circ}\text{C}$
EEO	: 0.15 - 0.2

### Risk and Safety Information

---

Storage Temperature( $^{\circ}\text{C}$ )	: Stored below $30^{\circ}\text{C}$
---	-------------------------------------

### Transport Information

---

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods