

RM9586 Gadolinium foil (Thickness 1mm, 25mm x 25mm)

Product Number Packing

RM9586 : 1NO

Product Information

Product Code : RM9586
Product Name : Gadolinium foil (Thickness 1mm, 25mm x 25mm)
Atomic Formula : G d
Atomic Weight : 157.25
CAS No. : 7440-54-2
EC No. : 231-162-2
HS Code : 2805 30 90
Shelf Life : 4 years

Technical Specification

Appearance : Grey silvery foil
Solubility : 33.3 mg soluble in 1 mL of dilute acids
Assay (EDTA Titration) : min. 99.00 %

GHS Safety Information

Hazard Statement(s) : H261
Precautionary Statement(s) : P231+P232-P422
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

WGK : 3
Storage Temperature(°C) : store below 30°C

Transport Information

UN No. : 3208
Class : 4.3
Packaging Group : 3
Marine Pollutant : No
ADR/RID : 3208 4.3/PG 3
IMDG : 3208 4.3/PG 3
IATA : 3208 4.3/PG 3