

AS031 Benzene, Hi-LR™

Product Number Packing

AS031	: 2.5L
AS031	: 500ML

Product Information

Product Code	: AS031
Product Name	: Benzene, Hi-LR™
Molecular Formula	: C ₆ H ₆
Molecular Weight	: 78.11
CAS No.	: 71-43-2
EC No.	: 200-753-7
HS Code	: 2902 20 00
EC Index No.	: 601-020-00-8
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear, Colourless Liquid.
Density (d 20/4)	: 0.878 - 0.879 g/mL
Water (K.F.)	: <= 0.1%
Non-volatile substances	: <= 0.002%
Assay (GC)	: min. 99.00 %

GHS Safety Information

Hazard Statement(s)	: H225-H304-H315-H319-H340-H350-H372
Precautionary Statement(s)	: P201-P210-P301+P310-P305+P351+P338-P308+P313-P331
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 45-46-11-36/38-48/23/24/25-65
S-Phrase(s)	: 53-45
WGK	: 3
RTECS	: CY1400000
Flash Point(°F)	: 12.2 °F
Flash Point(°C)	: -11 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information

UN No.	: 1114
Class	: 3



Product Information

Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	1114 3/PG 2
IMDG	:	1114 3/PG 2
IATA	:	1114 3/PG 2

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein