

1-Butanol, Hi-LR™

AS035


Product Identifier

CAS No.	:	71-36-3
EC No.	:	200-751-6
Molecular Formula	:	C ₄ H ₁₀ O
Molecular Weight	:	74.12
Synonym	:	n-Butyl alcohol
HS Code	:	2905 13 00
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	Clear colourless liquid
Identity (IR)	:	Passes test
Density (d 20/4)	:	~0.81 g/mL
Water (K.F.)	:	<= 0.20%
Non-volatile substances	:	<= 0.005%
Assay (GC)	:	min.99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H226- H302- H315- H318- H335- H336
Precautionary Statement(s)	:	P210- P280- P301+P312+P330- P304+P340+P312- P305+P351+P338+P310
UN No.	:	1120
Class	:	3
Packing Group	:	III
RTECS	:	EO1400000