

AS038 Chlorobezene, Hi-AR™/ACS

Product Number	Packing
----------------	---------

AS038	500ML
AS038	2.5L

Product Information

Product Code	: AS038
Product Name	: Chlorobezene, Hi-AR™/ACS
Synonym	: Monochlorobenzene; Benzene chloride
Molecular Formula	: C ₆ H ₅ Cl
Molecular Weight	: 112.56
CAS No.	: 108-90-7
EC No.	: 203-628-5
HS Code	: 2903 91 00
EC Index No.	: 602-033-00-1
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colorless liquid
Density (d 25/4)	: ~1.106 g/mL
Copper (Cu)	: <=0.00005%
Iron (Fe)	: <=0.00005%
Lead (Pb)	: <=0.00005%
Water (K.F.)	: <=0.05%
Non-volatile substances	: <=0.005%
Assay (GC)	: min. 99.50%

GHS Safety Information

Hazard Statement(s)	: H226-H332-H411
Precautionary Statement(s)	: P273
Signal Word	: Warning
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 10-20-51/53
S-Phrase(s)	: 24/25-61
WGK	: 2
RTECS	: CZ0175000
Flash Point(°F)	: 80.6 °F
Flash Point(°C)	: 27°C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	





Product Information

Transport Information

UN No.	:	1134
Class	:	3
Packaging Group	:	3
Marine Pollutant	:	YES
ADR/RID	:	1134 3/PG 3
IMDG	:	1134 3/PG 3
IATA	:	1134 3/PG 3

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein