

**AS041 Chloroform specially dried, Hi-AR™****Product Number Packing**

AS041 : 500ML

**Product Information**

Product Code : AS041  
Product Name : Chloroform specially dried, Hi-AR™  
Synonym : Trichloromethane; Methylidene trichloride  
Molecular Formula : CHCl<sub>3</sub>  
Molecular Weight : 119.38  
CAS No. : 67-66-3  
EC No. : 200-663-8  
HS Code : 2903 13 00  
Shelf Life : 3 years

**Technical Specification**

Appearance : Clear colourless liquid  
Solubility : 1 mL miscible with 1 mL of alcohol  
FTIR (Liquid film) : Matches with the standard pattern.  
Density (d 25/4) : 1.487 - 1.489 g/mL  
Water (K.F.) : <= 0.015%  
Free acid (as HCl) : <= 0.004%  
Ethanol (GC) : 0.6 - 1.0%  
Assay (GC) : 99.00 - 99.40%

**GHS Safety Information**

Hazard Statement(s) : H302+H332-H315-H319-H336-H351-H361d-H373  
Precautionary Statement(s) : P261-P281-P305+P351+P338  
Signal Word : Warning  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22-38-40-48/20/22  
S-Phrase(s) : 36/37  
WGK : 3  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

UN No. : 1888  
Class : 6.1  
Packaging Group : 3



## Product Information

Marine Pollutant	:	No
ADR/RID	:	1888 6.1/PG 3
IMDG	:	1888 6.1/PG 3
IATA	:	1888 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein