

## AS063 Petroleum ether, 40 - 60°C, Hi-AR™

### Product Number Packing

AS063	: 500ML
AS063	: 2.5L

### Product Information

Product Code	: AS063
Product Name	: Petroleum ether, 40 - 60°C, Hi-AR™
CAS No.	: 8032-32-4
EC No.	: 232-453-7
HS Code	: 2710 12 25
Shelf Life	: 3 years

### Technical Specification

Appearance	: Clear colourless liquid.
Density (d 20/4)	: ~0.65 g/mL
Water (K.F.)	: <=0.015%
Acidity	: <=0.0013 mL N%
Non volatile matter	: <=0.001%
Sulphur compound (CS <sub>2</sub> )	: <=0.0003%
Copper (Cu)	: <=0.0001%
Iron (Fe)	: <=0.0001%
Lead(Pb)	: <=0.0001%

### GHS Safety Information

Hazard Statement(s)	: H225-H304-H315-H336-H361-H373-H411
Precautionary Statement(s)	: P210-P261-P273-P281-P301+P310-P331
Signal Word	: Danger
Hazard Pictogram(s)	



### Risk and Safety Information

R-Phrase(s)	: 11-51/53-65-67
S-Phrase(s)	: 61-62
WGK	: 2
Flash Point(°F)	: -40 °F
Flash Point(°C)	: -40 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



### Transport Information

UN No.	: 1268
Class	: 3



## Product Information

Packaging Group	:	2
Marine Pollutant	:	Yes
ADR/RID	:	1268 3/PG 2
IMDG	:	1268 3/PG 2
IATA	:	1268 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : [info@himedialabs.com](mailto:info@himedialabs.com). Visit us at our website: [www.himedialabs.com](http://www.himedialabs.com)

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein