

**AS099 Acetonitrile, for pesticide residual analysis****Product Number Packing**

AS099 2.5L

**Product Information**

Product Code : AS099-2.5L  
Product Name : Acetonitrile, for pesticide residual analysis  
Synonym : Methyl cyanide, ACN  
Molecular Formula : C<sub>2</sub> H<sub>3</sub> N  
Molecular Weight : 41.05  
CAS No. : 75-05-8  
EC No. : 200-835-2  
HS Code : 2926 90 95  
EC Index No. : 608-001-00-3

**Technical Specification**

Appearance : Clear colourless liquid  
Density (d 25/4) : ~0.786 g/mL  
Water (K.F.) : <= 0.05%  
Pesticide analysis : <= 10 ng/l  
(Ethylparathion/PND)  
Pesticide analysis (Lindane/ECD) : <= 5 ng/l  
Assay (GC, on anhydrous substance) : min. 99.70 %

**GHS Safety Information**

Hazard Statement(s) : H225-H302-H312-H319-H332  
Precautionary Statement(s) : P210-P280-P305+P351+P338  
Signal Word : Danger  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 11-20/21/22-36  
S-Phrase(s) : 16-36/37  
WGK : 2  
RTECS : AL7700000  
Flash Point(°F) : 35.6 °F  
Flash Point(°C) : 2 °C  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**



## Product Information

UN No.	:	1648
Class	:	3
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	1648 3/PG 2
IMDG	:	1648 3/PG 2
IATA	:	1648 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein