

**AS103 Benzalkonium chloride (50 wt. % aqueous solution)****Product Number Packing**

AS103	: 5L
AS103	: 1L
AS103	: 500ML


**Product Information**

Product Code	: AS103
Product Name	: Benzalkonium chloride (50 wt. % aqueous solution)
Synonym	: Alkylbenzyltrimethylammonium chloride
HS Code	: 3402 12 00
Shelf Life	: 3 years


**Technical Specification**

Appearance	: Colourless to pale yellow liquid
Solubility	: 1 mL miscible in 1 mL of water
Assay (AT)	: 48.00 - 55.00%

**GHS Safety Information**

Hazard Statement(s)	: H302-H314-H400
Precautionary Statement(s)	: P273-P280-P305+P351+P338-P310
Signal Word	: Danger
Hazard Pictogram(s)	

**Risk and Safety Information**

R-Phrase(s)	: 21/22-34-50
S-Phrase(s)	: 26-36/37/39-45-61
WGK	: 2
RTECS	: BO3151000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

**Transport Information**

UN No.	: 3265
Class	: 8
Packaging Group	: 3
Marine Pollutant	: No
ADR/RID	: 3265 8/PG 3
IMDG	: 3265 8/PG 3
IATA	: 3265 8/PG 3