

2- Propanol, 69 - 71%

AS125

Product Identifier

CAS No.	:	67-63-0
EC No.	:	200-661-7
Molecular Formula	:	C ₃ H ₈ O
Molecular Weight	:	60.10
HS Code	:	2905 12 20
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	Clear colourless liquid
Solubility	:	1 mL miscible in 1 mL of water
Density (d 25/4)	:	0.853 - 0.863 g/mL
Water (K.F.)	:	25.0 - 30.0%
Concentration (GC)	:	69.00 - 71.00%

Safety Information 25.0 - 30.0%

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H336- H319- H225
Precautionary Statement(s)	:	P261- P280- P305+P351+P338- P210- P240
UN No.	:	1219
Class	:	3
Packing Group	:	II
RTECS	:	NT8050000