

## CG258 20 mL Makarthy Tube, Flat Bottom

### Product Number

### Packing

CG258 : 1 X 100 No.

### Product Information

Product Code : CG258  
Product Name : 20 mL Makarthy Tube, Flat Bottom

### Technical Specification

Description : Clear, tubular, flat bottom makarthy glass tube having screw neck made of borosilicate glass, Aluminium cap with silicone rubber gasket. The entire surface is smooth without presence of any air bubble.

Height of Tube : 84.0 - 86.0 mm

Body Diameter of Tube : 27.50 - 28.50 mm

Inner mouth Diameter of Tube : 21.80 - 22.20 mm

Outer mouth Diameter of Tube : 27.50 - 28.00 mm

Heat resistance test : On exposure to steam sterilization, Tube can withstand heat at 121°C for 30 minutes.

Weight with Cap : 23.00 - 25.00 gm

Overflow capacity : 35.0 - 40.0 mL

### Transport Information

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods