

## Test Tube Racks 'Z' Shape

CG259, CG260, CG261 & CG262

'Z' shape test tube racks are made from stainless steel with holes arranged in 5X5 and 5X10 format for different size tubes..

**Application :** For laboratory use as test tube racks.

Product Name	Product Code	Description
Test-Tube Racks 'Z' Shape	CG259	25-place test tube rack for 12-13mm diameter tubes.
	CG260	50-place test tube rack for 12-13mm diameter tubes.
	CG261	25-place test tube rack for 15-16mm diameter tubes.
	CG262	50-place .Test tube rack for 15-16mm diameter tubes.

### Product features :

- Stainless steel "Z" shape.
- Racks are rust proof and will hold up to repeated use.
- When not in use, the racks can be stacked inside one another.
- Pack size : 5X1 Nos.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.