



# Product Information

## CG281 0.5 mL Autoclavable PCR Tubes, Flat Lid

### Product Number

### Packing

CG281 : 1 x 1000 No.

### Product Information

Product Code : CG281

Product Name : 0.5 mL Autoclavable PCR Tubes, Flat Lid

### Technical Specification

Description : Autoclavable PCR tubes having conical base with flat lid and with graduation. It is made up of polypropylene material.

Clarity : Clear without any particles and fibers.

Height of Tube : 30.0 - 31.0. mm

Height of Cap : 4.0 - 5.0 mm

Inner mouth Diameter of tube : 6.5 - 7.5 mm

Outer mouth Diameter of tube : 9.0 - 9.5 mm

Outer mouth Diameter of cap : 6.5 - 7.0 mm

Capacity : 0.5 mL

Temperature tolerance test : Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.

Product alignment : Lid closer of the tube is tight.

### Transport Information

ADR/RID : Not Dangerous Goods

IMDG : Not Dangerous Goods

IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein