

CG282 0.2 mL Autoclavable PCR Tubes, Flat Lid

Product Number	Packing
----------------	---------

CG282	: 1 x 1000 No.
-------	----------------

Product Information

Product Code	: CG282
Product Name	: 0.2ml PCR Tube, PP material, Conical Bott

Technical Specification

Description	: Autoclavable PCR tubes having conical base with flat lid. It is made up of polypropylene material.
Clarity	: Clear without any particles and fibers.
Height of Tube	: 19.5 - 20.5 mm
Height of Cap	: 3.5 - 4.5 mm
Inner mouth Diameter of tube	: 5.0 - 6.0 mm
Outer mouth Diameter of tube	: 6.5 - 7.5 mm
Outer mouth Diameter of cap	: 5.0 - 6.0 mm
Capacity	: 0.2 mL
Temperature tolerance test	: Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Product alignment	: Lid closer of the tube is tight.

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods