

CG283 0.2 mL Autoclavable PCR Tubes, Dome shaped Lid

Product Number	Packing
----------------	---------

CG283	: 1 x 1000 No.
-------	----------------

Product Information

Product Code	: CG283
Product Name	: 0.2 mL AutoclavablePCRTubes, Dome shaped Lid

Technical Specification

Description	: Autoclavable PCR tubes having conical base with dome shaped lid. It is made up of polypropylene material.
Clarity	: Clear without any particles and fibers.
Height of Tube	: 20.5 - 21.5 mm
Height of Cap	: 5.0 - 6.0 mm
Inner mouth Diameter of tube	: 5.0 - 6.0 mm
Outer mouth Diameter of tube	: 6.5 - 7.5 mm
Outer mouth Diameter of cap	: 5.0 - 6.0 mm
Capacity	: 0.2 mL
Temperature tolerance test	: Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Product alignment	: Lid closer of the tube is tight.

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.