



Product Information

CG300

96-places (8x12) Autoclavable Orange colour PCR Tube Workstation

Product Number**Packing**

CG300 : 1 x 5 No.

Product Information

Product Code : CG300
Product Name : 96-places (8x12) Autoclavable Orange colour PCR Tube Workstation

Technical Specification

Description : Autoclavable Orange colour PCR tube workstation with alphanumerically grid marking, it is made up of polypropylene material.
Dimension : 215 x 112 x 27 mm (L x W x H)
Clarity : Clear without any scratches and spots.
Product alignment : Proper grid marking of alphanumerical figure.
Capacity : 96 Place (8x12) for 0.2 ml, 0.5 ml & 1.5 ml PCR tubes
Temperature tolerance test : Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein