

CG305

1000 µL Autoclavable Blue Colour Pipet Tips

Product Number

Packing

CG305 : 1 x 500 No.

Product Information

Product Code : CG305
Product Name : 1000 µL Autoclavable Blue Colour Pipet Tips

Technical Specification

Description : Autoclavable yellow coloured pipet tips made up of polypropylene material.
Clarity : Clear transparent without any particles and Fibers.
Inner Dimension : Approx 7.5 mm
Capacity : 1000 µl
Product alignment : Straight without any bending.
Temperature tolerance test : Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Compatible Pipette Code : LA246, LA247, LA248, LA606, LA607, LA608, LA614, LA616.
Use : Used to dispense liquid samples.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods