

CMS1673**Dapson****Product Number****Packing**

CMS1673 : 5G

Product Information

Product Code : CMS1673
Product Name : Dapson
Synonym : 4,4'-Diaminodiphenyl sulphone; 4-Aminophenyl sulphone
Molecular Formula : C₁₂H₁₂N₂O₂S
Molecular Weight : 248.30
CAS No. : 80-08-0
EC No. : 201-248-4
HS Code : 2930 90 99
EC Index No. : 612-084-00-1
Shelf Life : 3 years

Technical Specification


Appearance : White to tan crystals or powder
Solubility : 33.3 mg soluble in 1 mL of methanol
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 175 - 181°C
Assay (HPLC/NT, on dry basis) : 97.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H302
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22
S-Phrase(s) : 22
WGK : 1
RTECS : BY8925000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods



Product Information

IATA

: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein