

CMS4567 4-Aminosalicylic acid

Product Number Packing

CMS4567 : 5G

Product Information

Product Code : CMS4567
Product Name : 4-Aminosalicylic acid
Molecular Formula : $C_7H_7NO_3$
Molecular Weight : 153.14
CAS No. : 65-49-6
EC No. : 200-613-5
HS Code : 2922 50 00
Shelf Life : 3 years

Technical Specification

Appearance : White to grey or brown powder
Solubility : 33.3 mg soluble in 1 mL of ethanol
FTIR : Matches with the standard pattern
Melting range : 130 - 157°C
Assay (NaOH Titration/NT) : 98.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H319
Precautionary Statement(s) : P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 22-36
S-Phrase(s) : 26
WGK : 2
RTECS : VO1225000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods