

CMS527 Demeclocycline hydrochloride**Product Number Packing**

CMS527 : 1G

Product Information

Product Code : CMS527
Product Name : Demeclocycline hydrochloride
Synonym : 7-Chloro-6-demethyltetracycline hydrochloride
Molecular Formula : $C_{21}H_{21}ClN_2O_8 \cdot HCl$
Molecular Weight : 501.31
CAS No. : 64-73-3
EC No. : 200-592-2
HS Code : 2941 30 00
Shelf Life : 3 years

Technical Specification

Appearance : Yellow powder or solid
Solubility : 10 mg soluble in 1 mL of water
Assay (HPLC) : min. 99.00%

GHS Safety Information

Hazard Statement(s) : H317
Precautionary Statement(s) : P280
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 43
S-Phrase(s) : 36/37
WGK : 2
RTECS : QI7700000
Storage Temperature(°C) : Store at 2 - 8°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods