

EP6 96 well, Medium binding, detachable, 8x12 Elisa plates

Product Number	Packing
EP6	: 5 x 10 No.
EP6	: 10 x 10 No.
EP6	: 20 x 10 No.

Product Information

Product Code	: EP6
Product Name	: 96 well, Medium binding, detachable, 8x12 Elisa plates

Technical Specification

Description	: Medium binding capacity of detachable Elisa plates having alphanumerically grid marking, it is made up of polystyrene material.
Capacity	: 96 well (8x12)
Clarity	: Clear translucent without any scratches, spots and breakage.
Product alignment	: Proper grid marking of alphanumeric figure.

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.