



Technical Data

Campylobacter Growth Supplement

FD009

A chemical enrichment recommended for the enhancement of growth and aerotolerance of *Campylobacter fetus* .

Composition

Per vial sufficient for 500 ml medium

Ingredients

Sodium pyruvate
Sodium metabisulphite
Ferrous sulphate

Concentration

0.125g
0.125g
0.125g

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water and mix well to dissolve. Avoid frothing of the solution. Aseptically add the contents to 500 ml of sterile, molten, cooled (45-50°C) Columbia Blood Agar Base [M144](#) / [M144A](#) / Columbia Blood Agar Base, HiVeg™ [MV144](#) / [MV144A](#) / Blood Agar Base No. 2 [M834](#) / [M834A](#) Blood Agar Base No. 2, HiVeg™ [MV834](#) / [MV834A](#) / Brucella Agar Base [M074](#) / [M074A](#) / Brucella HiVeg™ Agar Base [MV074](#) / [MV074A](#) / Brucella Broth Base [M348](#) / Brucella HiVeg™ Broth Base [MV348](#) with 5-7% v/v lysed defibrinated horse or sheep blood along with Campylobacter Supplement I [FD006](#) or Campylobacter Supplement II [FD007](#) or Campylobacter Supplement III [FD008](#) . Mix gently and pour into sterile petri plates

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Revision : 1 / 2012

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.