



# Technical Data

## Vancko T Supplement

FD028

An antibiotic supplement recommended for the selective isolation of *Neisseria gonorrhoeae* and *Neisseria meningitidis*.

### Composition

Per vial sufficient for 500 ml medium

| *Ingredients         | Concentration |
|----------------------|---------------|
| Vancomycin           | 1mg           |
| Colistin sulphate    | 3.750mg       |
| Amphotericin B       | 0.500mg       |
| Trimethoprim lactate | 1.500mg       |

### Directions:

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water. Mix gently. Avoid frothing of the solution to dissolve. Aseptically add the contents of one vial to 490 ml sterile, molten, cooled (45-50°C) GC Agar Base [M434](#) / GC HiVeg™ Agar Base [MV434](#) along with Haemoglobin Solution [FD022](#) and GC Supplement w/ Antibiotics [FD021](#). Also add sterile rehydrated contents of 1 vial of Yeast Autolysate Supplement [FD027](#) or Vitamino Growth Supplement [FD025](#). Mix well and pour into sterile petri plates.

### Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

\* Not For Medicinal Use

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.