



# Technical Data

## Legionella Selective Supplement IV (MWY)

FD040

An antibiotic supplement recommended for the selective isolation of *Legionella* species.

### Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Polymyxin B sulphate	25000Unit
Glycine	1.500g
Anisomycin	40mg
Vancomycin	0.500mg
Bromo thymol blue	5mg
Bromo cresol purple	5mg

### Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) Buffered Charcoal Yeast Extract Agar Base [M813](#) / Modified Buffered Charcoal Agar Base [M892](#) / Modified Buffered Charcoal HiVeg™ Agar Base [MV892](#) / Buffered Charcoal Yeast Extract Agar Base, Modified [M813I](#) alongwith rehydrated contents of 1 vial of Legionella Supplement(Twin pack) [FD041A](#) . Mix well and pour into sterile petri plates.

### Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Revision : 03 / 2019

\* Not For Medicinal Use

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.