



Technical Data

Mycoplasma Enrichment Supplement

FD075

A selective supplement for the isolation of Mycoplasma.

Composition

Per vial sufficient for 70 ml medium

*Ingredients	Concentration
Horse serum	20ml
Yeast extract (25% w/v solution)	10ml
Thallos acetate	25mg
Penicillin G	20000IU

Directions:

Aseptically add the contents of 1 vial to 70 ml of sterile, molten, cooled (45-50°C) Mycoplasma Agar Base [M266](#) / Mycoplasma HiVeg™ Agar Base [MV266](#) / Mycoplasma Broth Base w/o CV [M267](#) / Mycoplasma HiVeg™ Broth Base w/o CV [MV267](#) / Mycoplasma Broth Base w/ CV [M268](#) / Mycoplasma HiVeg™ Broth Base w/ CV [MV268](#) / PPLO Modified Broth Base w/o CV [M1586](#) . For Mycoplasma Broth Base w/ CV [M268](#) / Mycoplasma HiVeg™ Broth Base w/ CV [MV268](#) add 0.285 ml of 1% Potassium Tellurite [FD052](#) along with Mycoplasma Enrichment Supplement [FD075](#) . Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.