



Technical Data

Trichomonas Selective Supplement I

FD099

An antimicrobial supplement recommended for selective isolation of *Trichomonas* species from clinical specimens.

Composition

Per vial sufficient for 500 ml medium

*Ingredients

Streptomycin

Penicillin G

Concentration

500mg

125000Unit

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and avoid frothing of solution. Aseptically add the contents to 475 ml of sterile, molten, cooled (45-50°C) Kupferberg Trichomonas Broth Base (Trichomonas Broth Base, Kupferberg) [M305](#) / Kupferberg Trichomonas HiVeg™ Broth Base [MV305](#) along with 25 ml sterile beef or human serum. Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.