



Technical Data

Enterococcus Selective Supplement

FD148

An antibiotic supplement recommended for enumeration of intestinal Enterococci from surface and waste water.

Composition

Per vial sufficient for 100 ml medium

*Ingredients

Nalidixic acid
Thallos acetate

Concentration

25mg
200mg

Directions:

Rehydrate the contents of one vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to sterile, cooled (45-50°C) 100 ml of MUD SF Broth Base [M1343](#) along with 1 ml of sterile 1% T.T.C. Solution [FD057](#). Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.