



Technical Data

BSK-H Supplement

FD180

A sterile filtered rabbit serum to support growth of *Borrelia burgdorferi* .

Ingredients

Rabbit serum

Concentration

100ml

Directions:

To prepare complete medium warm up the refrigerated rabbit serum to room temperature and aseptically add 30 ml to 500 ml sterile, cooled (45-50°C) BSK-H Medium Base [M1668](#) / [LQ001](#) / BSK-H Medium Powder [M5107](#) along with one vial of Antibiotic Mixture of Borrelia [FD179](#) . Mix well and dispense as desired.

Storage and Shelf Life

Store at -20°C. Use before the expiry date on the label.

Revision : 3 / 2018

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.