



Technical Data

Albumin Glucose Supplement

FD201

An enrichment supplement recommended for rapid cultivation of *Mycobacterium tuberculosis*

Composition

Per vial sufficient for 200 ml medium

Ingredients

Albumin fraction V
Dextrose
Distilled water

Concentration

1g
1.500g
20ml

Directions:

Warm up the refrigerated supplement to room temperature. Aseptically add to sterile, cooled (45-50°C) 180 ml Dubos Broth Base [M067](#) / Dubos HiVeg™ Broth Base [MV067](#) . Mix well and pour into sterile screw capped tubes.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.