



# Technical Data

## Vancomycin Supplement

FD261

This antibiotic supplement is recommended for selective isolation of Vancomycin Resistant Enterococci.

### Composition

Per vial sufficient for 500 ml medium

#### \*Ingredients

Vancomycin

#### Concentration

3mg

### Directions:

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add 2 vials to 1000 ml of sterile molten Vancomycin Resistant Enterococci (VRE) Agar Base [M1763](#) along with 1 vial of Meropenem Supplement [FD262](#) . Mix well and pour into sterile Petri plates.

### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision : 1 / 2012

\* Not For Medicinal Use

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.