

GRM126 Cellulose powdered

Product Number Packing

GRM126 : 500G

Product Information

Product Code : GRM126
Product Name : Cellulose powdered
Synonym : α -Cellulose; Cotton linters
CAS No. : 9004-34-6
EC No. : 232-674-9
HS Code : 3912 90 90
Shelf Life : 4 years

Technical Specification

Appearance : White to off-white powder
Solubility : Insoluble in water
pH (10% aqueous slurry at 25°C) : 5.0 - 7.5
Identification Test : Passes
Loss on drying (at 105°C; 2hr) : \leq 7.0%
Chloride (Cl) : \leq 0.05%
Ash : \leq 0.3%
Water soluble substances : \leq 1.5%
Assay (FAS Titration, on dry basis) : 97.00 - 102.00%

Risk and Safety Information

WGK : 3
RTECS : FJ5691460
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods