

**GRM2033**

**Carbopol - 940**

**Product Number**

**Packing**

GRM2033 : 1KG  
GRM2033 : 500G

**Product Information**

Product Code : GRM2033  
Product Name : Carbopol - 940  
HS Code : 3906 90 90  
Shelf Life : 4 years

**Technical Specification**

Appearance : White hygroscopic powder  
Solubility : 1% w/v aqueous solution gives a slight opalescent gel at room temperature.  
Loss on drying (at 80°C, 2 hr) : <= 3.00%  
Viscosity : 40000 - 60000 cPs (0.5% neutral solution, spindle no. 7, 20 rpm at 25°C)  
Assay (Carboxylic acid, NaOH Titration) : 56.00 - 68.00%

**Risk and Safety Information**

Storage Temperature(°C) : Store below 30°C

**Transport Information**

Marine Pollutant : -  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods