



Product Information

RM7402 Polyethylene glycol 8000

Product Number **Packing**

RM7402 : 500G
RM7402 : 5KG

Product Information

Product Code : RM7402
Product Name : Polyethylene glycol 8000
Molecular Weight : ~8000
CAS No. : 25322-68-3
Other Information : For synthesis
Shelf Life : 4 years

Technical Specification

Appearance : White to creamy white flakes.
Solubility : Soluble in water.
pH (5% in water at 25°C) : 4.00 - 7.50
Sulphated ash : <= 0.00001%
Arsenic (As) : <= 0.0003%
Heavy metals : <= 0.0005%
Viscosity (@98.9°C) : 470-900 CST
Congealing point : 58-62°C

Risk and Safety Information

Storage Temperature(°C) : Store below 30°C