

**GRM7495 Sodium bromate, Hi-LR™**

**Product Number Packing**

---

GRM7495 : 500G

**Product Information**

---

Product Code : GRM7495  
Product Name : Sodium bromate, Hi-LR™  
Molecular Formula :  $\text{NaBrO}_3$   
Molecular Weight : 150.89  
CAS No. : 7789-38-0  
EC No. : 232-160-4  
HS Code : 2829 90 40  
Shelf Life : 4 years

**Technical Specification**

---

Appearance : White crystals or powder or chunks  
Solubility : 100 mg soluble in 1 mL of water  
Sulphate (SO<sub>4</sub>) : ≤ 0.02%  
Assay (Iodometry) : 99.00 - 101.00 %

**GHS Safety Information**

---

Hazard Statement(s) : H272-H302-H315-H319-H335  
Precautionary Statement(s) : P220-P261-P305+P351+P338  
Signal Word : Danger  
Hazard Pictogram(s)



**Risk and Safety Information**

---

R-Phrase(s) : 8-22-36/37/38  
S-Phrase(s) : 17-26-36  
WGK : 3  
RTECS : EF8750000  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)



**Transport Information**

---

UN No. : 1494  
Class : 5.1  
Packaging Group : 2  
Marine Pollutant : No  
ADR/RID : 1494 5.1/PG 2



## Product Information

IMDG : 1494 5.1/PG 2  
IATA : 1494 5.1/PG 2

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.  
email:info@himedialabs.com. Visit us at our website: [www.himedialabs.com](http://www.himedialabs.com)

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein