

GRM8535 Jeffamine® T-403**Product Number Packing**

GRM8535 : 500G

Product Information

Product Code : GRM8535
Product Name : Jeffamine® T-403
Molecular Formula : $C_2H_5C[CH_2[OCH_2CH(CH_3)]_nNH_2]_3$
Molecular Weight : ~ 440
CAS No. : 39423-51-3
HS Code : 3907 20 99
Other Information : Key: Jeffamine is a registered trademark of Huntsman Corporation
Shelf Life : 4 years

Technical Specification

Appearance : Colourless to pale yellow with slight haze liquid
Solubility : 1 mL miscible in 1 mL of water
pH (5% in water at 25°C) : 11.00 - 12.00
Primary amine, % of total amine (NT) : min. 90.00%

GHS Safety Information

Hazard Statement(s) : H302-H312-H314
Precautionary Statement(s) : P280-P305+P351+P338-P310
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 20/21/22-34
S-Phrase(s) : 23-26-27-36/37/39-45
WGK : 3
Flash Point(°F) : 235.4 °F
Flash Point(°C) : 113 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 3267
Class : 8
Packaging Group : 3
Marine Pollutant : No
ADR/RID : 3267 8/PG 3



Product Information

IMDG : 3267 8/PG 3
IATA : 3267 8/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein