

## Disposable Micro Tips for Autopipettes

**LA270, LA271, LA272,  
LA273, LA282, LA284,  
LA286, LA288, LA887  
& LA888.**

Made of polypropylene. Specifically designed to suit Autopipettes. These can be used on any regular micro pipette. Available in 3 different sizes: 1.00 to 10000.00 µl with different Color .

**Application :** Pipetting Applications, precious application, Molecular biology application including DNA sequencing and plasmid purification, PCR preparation, automated sample preparation and for use in a wide variety of pipetting applications etc.

Product Name	Product Code	Description	Volume Capacity (µl)	Materials
Disposable micro Tip-P.P.	LA270	Natural Color micro-tips.	10000.00	Polypropylene
Disposable micro Tip-P.P.	LA271	Natural Color suitable for Fixvo and variable micropipettes, Compatible Pipette codes are LA609, LA610 & LA619.	5000.00	Polypropylene
Disposable micro Tip-P.P.	LA272	Natural Color suitable for Harmony, Fixvo and variable micropipette, Compatible Pipette codes are LA246, LA247, LA248, LA606, LA607, LA608, LA614, LA616, LA742, LA748 & LA776.	1000.00	Polypropylene
Disposable micro Tip-P.P.	LA273	Pink Color suitable for Nichiryu auto pipette / Fixvo micropipette , Compatible Pipette codes are LA629, LA630 & LA631.	1000.00	Polypropylene
Disposable micro Tip-P.P.	LA282	Green Color suitable for variable micropipette , Compatible Pipettes codes are LA629.	250.00	Polypropylene
Disposable micro Tip-P.P.	LA284	Yellow Color suitable for harmony and fixvo II micropipette, Compatible Pipettes codes are LA235, LA237, LA240, LA242, LA599, LA600, LA601, LA602, LA603, LA604, LA605, LA612, LA613, LA615, LA618, LA627, LA628, LA739, LA741, LA745, LA746, LA747, LA778, LA779, LA781, LA791, LA797 & LA798.	200.00	Polypropylene
Disposable micro Tip-P.P.	LA286	Natural Color suitable for Fixvo , harmony and variable micropipette, Compatible Pipettes codes are LA235, LA237, LA240, LA242, LA599, LA600, LA601, LA602, LA603, LA604, LA605, LA612, LA613, LA615, LA618, LA627, LA628, LA739, LA741, LA745, LA746, LA747, LA778, LA779, LA781, LA791, LA797 & LA798	200.00	Polypropylene
Disposable micro Tip-P.P.	LA288	Natural Color suitable for Clickpet Dispensers and Fixvo, Harmony and variable micropipettes. Compatible Pipette codes are LA233, LA235, LA237, LA240, LA623, LA624, LA611, LA617, LA738, LA743, LA744, LA775, LA777 & LA780.	1.00 to 20.00	Polypropylene
Disposable micro Tip-P.P.	LA887	Natural colour with long narrow tip, Compatible Pipette code is LA878.	10000	Polypropylene
Disposable micro Tip-P.P.	LA888	Natural colour with long narrow tip. Compatible Pipette code is LA878.	10000	Polypropylene

**Product Features:**

- Disposable tips made of polypropylene.
- Available in different sizes and colors.
- Pack Size: 500 Nos and 100 Nos . in a pack.

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

---

HiMedia Laboratories Pvt. Ltd. A-516,Swastik Disha Business Park,Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India.  
Customer care No.: 022-6147 1919 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)