



Product Information

LA321 Indicator pH paper Range 6.5 - 9.0

Product Number

Packing

LA310 : 1PK

Product Information

Product Code : LA321
Product Name : Indicator pH paper Range 6.5 - 9.0
Shelf Life : 1 year

Technical Specification

Description : Yellow coloured strip of filter paper booklet describing pH range 6.5 - 9.0
Result : At pH 6.5 colour of the filter paper strip remains yellow and at pH 9.0 colour changes to violet

Risk and Safety Information

Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : -
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods