

## e-Stand

LA554

e-Stand is pipette stand. e-Stand is blue in color with 10 spaces, and two to three sturdy feet in size. It can be stored vertically. It is lightweight, comfortable to use. Constructed of heavy, solvent-resistant plastic and holds a maximum of 10 Pipettes.

### Pipette Protection

- It offers strong resistance to aggressive chemicals and UV light
- Protects valuable instruments from shocks and contamination.
- Pipettes hang vertically to prevent liquids from running inside the shaft.

### Comfortable Operation

- Its rotating head gives an easy access to the pipette of choice and
- Adapted to limited work bench space

Product Name	Product Code	Description	Capacity
e-Stand	LA554	Pipette Stand	10 Pipettes

### Product Features

1. Autoclavable
2. Packaging: 1 No.

#### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.