



Product Information

LA612

Varivol II Micropipette - 50

Product Number

Packing

LA612 : 1 NO

Product Information

Product Code : LA612
Product Name : Varivol II Micropipette - 50

Technical Specification

Description : Yellow coloured Variable volume microlitre pipettes wherein volumes can be continuously adjusted by rotating the plunger or thumb button. It has click-stop digital system, reverse pipetting technique and separate, streamlined tip ejector. Each pipette has been tested, calibrated and comes with its own calibration certificate.

Capacity : 5 - 50 μ L

Use : Used in laboratory to dispense liquid samples.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein