

## Stainless Steel Forceps

**LA706, LA707, LA708,  
LA709, LA710, LA711,  
LA721, LA767.**

These Stainless Steel Forceps are straight with blunt and pointed tips and handles. Available in three different sizes.

These high quality and durable forceps have great picking up quality. These can be used for slide making, for handling items safely etc.

**Application :** Ideal for handling and moving membrane to and from filter holders and Petri dishes. Stainless-steel forceps are also used in plant tissue culture work, grasping of tissues or dressings etc.

Product Name	Product Code	Description	Size (Inch)
Stainless -Steel - Forceps	LA706	Blunt Stainless Steel Forceps	6.00
	LA707	Blunt Stainless Steel Forceps	8.00
	LA708	Blunt Stainless Steel Forceps	10.00
	LA709	Pointed Stainless Steel Forceps	6.00
	LA710	Pointed Stainless Steel Forceps	8.00
	LA711	Pointed Stainless Steel Forceps	10.00
	LA721	Blunt Stainless Steel Forceps	9.00
	LA767	Blunt Stainless Steel Forceps	12.00

### Product Features:

- High-quality, durable stainless steel
- Easy to handle.
- Autoclavable.
- Pack Size : 2 Nos in a pack

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.