



# Product Information

**LA749**

**Barrier Tips, Capacity: 10 µl**

---

**Product Number****Packing**

---

LA749 : 1 x 96 No.  
LA749 : 10 x 96 No.

---

**Product Information**

---

Product Code : LA749  
Product Name : Barrier Tips, Capacity: 10 µl

---

**Technical Specification**

---

Description : Neutral colour pre-sterile universal low retention barrier tips presence of filter. Thin wall design provides an optimal seal for single and multichannel pipettors. Tips are Nuclease and endotoxin free.

Clarity : Clear transparent without any particles and Fibers.

Inner Dimension : Approx 3.5 mm

Capacity : 10 µl

Product alignment : Straight without any bending.

Compatible Pipette Code : LA233, LA235, LA237, LA240, LA623, LA624, LA611, LA617.

Use : Used to dispense small amounts of liquid.

---

**Transport Information**

---

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein