



Product Information

LA751

Barrier Tips, Capacity: 200 µl

Product Number**Packing**

LA751 : 1 x 96 No.
LA751 : 10 x 96 No.

Product Information

Product Code : LA751
Product Name : Barrier Tips, Capacity: 200 µl

Technical Specification

Description : Neutral colour pre-sterile universal low retention barrier tips presence of filter. Thin wall design provides an optimal seal for single and multichannel pipettors. Tips are Nuclease and endotoxin free.

Clarity : Clear transparent without any particles and Fibers.

Inner Dimension : Approx 5.3 mm

Capacity : 200 µl

Product alignment : Straight without any bending.

Compatible Pipette Code : LA235, LA237, LA240, LA242, LA599, LA600, LA601, LA602, LA603, LA604, LA605, LA612, LA613, LA615, LA618, LA627, LA628.

Use : Used to dispense small amounts of liquid.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods