



Product Information

LA859

Barrier Tips, Capacity: 1000 µl

Product Number

Packing

LA859	: 1 x 96 No.
LA859	: 8 x 96 No.

Product Information

Product Code	: LA859
Product Name	: Barrier Tips, Capacity: 1000 µl

Technical Specification

Description	: Neutral colour pre-sterile universal low retention barrier tips presence of filter. Thin wall design provides an optimal seal for single and multichannel pipettors. Tips are Nuclease and endotoxin free.
Clarity	: Clear transparent without any particles and Fibers.
Inner Dimension	: Approx 7.5 mm
Capacity	: 1000 µl
Product alignment	: Straight without any bending.
Compatible Pipette Code	: LA246, LA247, LA248, LA606, LA607, LA608, LA614, LA616.
Use	: Used to dispense liquid samples.

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein