



Product Information

LA886 Bluple Nitrile Examination Gloves, Medium Size

Product Number

Packing

LA886	: 10 X 50 No.
LA886	: 10 X 100 No.

Product Information

Product Code	: LA886
Product Name	: Bluple Nitrile Examination Gloves, Medium Size

Technical Specification

Description	: Bluple medium sized synthetic nitrile powder free examination gloves having finger texture and fully textured surface for a better grip. The gloves are longer and thicker for greater protection over the wrist and arm.
Weight	: 3.0 - 4.0 gm

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein