



# Product Information

## LA894 Bluple Nitrile Examination Gloves, Large Size

### Product Number

### Packing

---

LA894	: 10 X 50 No.
LA894	: 10 X 100 No.

### Product Information

---

Product Code	: LA894
Product Name	: Bluple Nitrile Examination Gloves, Large Size

### Technical Specification

---

Description	: Bluple large sized synthetic nitrile powder free examination gloves having finger texture and fully textured surface for a better grip. The gloves are longer and thicker for greater protection over the wrist and arm.
Weight	: 3.0 - 4.0 gm

### Transport Information

---

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein