



# Product Information

**LA967       $\mu$ Pet Autoclavable Micropipette Capacity: 8 Channel, (30 - 300  $\mu$ L)**

---

<b>Product Number</b>	<b>Packing</b>
-----------------------	----------------

---

LA967	: 1 No.
-------	---------

---

### Product Information

---

Product Code	: LA967
Product Name	: $\mu$ Pet Autoclavable Micropipette Capacity: 8 Channel, (30 - 300 $\mu$ L)

---

### Technical Specification

---

Description	: Fully autoclavable $\mu$ pet multi channel micropipette wherein volume can be adjusted by using the variable volume controller knob which lock or unlock the volume dial by pushing the lock lever upwards and downwards. It has reverse pipetting technique and separate streamlined tip ejector. Each pipette has been tested and calibrated and comes with its own
Capacity	: 8 Channel, (30 - 300 $\mu$ L)
Use	: Used in laboratory to dispense liquid samples.

---

### Transport Information

---

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein