

LA972 200 µL Yellow colour Micropipette Tips

Product Number

Packing

LA972	: 1 X 1000 No.
LA972	: 5 X 1000 No.

Product Information

Product Code	: LA972
Product Name	: 200 µL Yellow colour Micropipette Tips

Technical Specification

Description	: Autoclavable Yellow coloured clear nuclease free, pyrogen free micropipette tips.
Clarity	: Clear transparent without any particles and Fibers.
Inner Dimension	: Approx 5.3 mm
Capacity	: 200 µL
Product alignment	: Straight without any bending.
Compatible Pipette Code	: LA599, LA600, LA604, LA612, LA613, LA615, LA618, LA959, LA966.
Temperature tolerance test	: Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Use	: Used to dispense small amount of liquid.

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods