

LA973 Micropipette Tips Capacity : 1000 µL

Product Number

Packing

LA973	: 1 X 1000 No.
LA973	: 5 X 1000 No.

Product Information

Product Code	: LA973
Product Name	: Micropipette Tips Capacity : 1000 µL

Technical Specification

Description	: Autoclavable Blue coloured clear nuclease free, pyrogen free micropipette tips.
Clarity	: Clear transparent without any particles and Fibers.
Inner Dimension	: Approx 7.4 mm
Capacity	: 1000 µL
Product alignment	: Straight without any bending.
Compatible Pipette Code	: LA606, LA607, LA608, LA614, LA616, LA960, LA963.
Temperature tolerance test	: Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Use	: Used to dispense small amount of liquid.

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods